



## Judd Rice Laboratory

at the USC/Norris Comprehensive Cancer Center

[www.histonecode.com](http://www.histonecode.com)

## Immunofluorescence

Revised by J Sims  
June 2004

### Procedure

---

1. Plate cells on coverslips to a density of 105 cells/well (2mL/well) in 6 well dishes.
2. Fix cells when they reach ~60-80% confluent.
  - 2.1. Make 4% paraformaldehyde/PBS- dissolve 4g paraformaldehyde in 100mL PBS, heat to 60°C for it to go into solution.
  - 2.2. Add 2mL paraformaldehyde/well and incubate 10 min at RT.
  - 2.3. Aspirate paraformaldehyde and wash 3x10min with PBS
3. Permeabilize cells with 0.2% Triton X-100/PBS (200µL Triton/100mL PBS).
4. Incubate 5 min RT.
5. Aspirate Triton and wash 3x5min with PBS
6. Block cells in 5% donkey serum/PBS (5mL donkey serum/100mL PBS) for 1 hour at RT.
7. Aspirate serum.
8. Add primary Ab, diluted in 5% donkey serum/PBS (125µL/coverslip) to the center of the coverslip. Try to move the coverslip to the center of the well to prevent the antibody from wicking off of the coverslip.
9. Incubate at 37°C for 1 hour in humidifier.  
(Or a Tupperware container lined with damp paper towels)
10. Add 2mL 5% donkey serum and wash 5 min (repeat 2 more times).
11. Aspirate serum and add secondary Ab, diluted in 5% donkey serum/PBS (125µL/coverslip) to the center of the coverslip.
12. Incubate at 37°C for 1 hour in humidifier.
13. Wash 3x10 min with 5% donkey serum/PBS.

14. Mount coverslips using mounting medium containing DAPI.
  - 14.1. Add one drop of mounting medium to slide.
  - 14.2. Rinse coverslip in dH<sub>2</sub>O.
  - 14.3. Place coverslip cells side down onto mounting medium.
  - 14.4. Blot off edges and let dry.
  
15. Store slides at 4°C